BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFIA/BIBLIOGRAFIA

- Panteghini, M. & Blais, R., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
- Dufour, D.R., *The Pancreas: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 651 and appendix.
- Winn-Deen, E.S., *et al. Clin. Chem.*, (1988), **34**, 2005
- Guder W.G, *et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma, and serum samples*, (2002). WHO/DIL/LAB/99.1 Rev.2
- Schumann, G., *et al. Clin. Chem. Lab Med.*, (2006), **44**, 1146.
- Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACCPress, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACCPress, (1995).

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
- Symbols used on our documentation are defined on ISO-15223-1 standard, except for some presented in the symbols glossary available on the VitalScientific Website. (Symbols glossary).

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- Special Programming Instructions :
see § INSTALLATION AND USE

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English - EN

INTENDED USE

AMYLASE SL is an *in vitro* diagnostic reagent intended for the quantitative determination of amylase in human serum and plasma samples on analyzers. This *in vitro* diagnostic device is for professional use only.

CLINICAL SIGNIFICANCE ⁽¹⁻²⁾

α -Amylase is an enzyme mainly from pancreatic or salivary origin that hydrolyses 1,4- α -glycosidic bonds, thus helping for starch digestion. An increase of serum amylase can be observed in acute pancreatitis mainly, as well as in many other situations such as other pancreatic or intra-abdominal pathologies, and salivary gland diseases. During acute pancreatitis, serum amylase level starts rising after a few hours and can be increased up to 6-fold, peaking after 12 to 72h, and returning to normal level within about 4 days. In clinical practice, measurement of amylase activity is mainly used to help diagnose and monitor acute pancreatitis, and may also help diagnose and follow-up other pancreatic diseases.


LIMITATION OF USE

The quantitative assay of amylase alone cannot be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

METHOD & PRINCIPLE ⁽³⁾

CNP-G₃ / Standardization against IFCC method – Kinetic

CNP-G₃ substrate is hydrolyzed by catalytic action of α -amylase to produce CNP (2-chloro-4-nitrophenol) which absorbs at 405 nm.

CNP-G₃ = 2-Chloro-4-nitrophenyl- α -maltotrioside
CNP-G₂ = 2-Chloro-4-nitrophenyl- α -maltoside

The rate of increase in absorbance is measured at 405 nm and is directly proportional to the activity α -amylase in the sample.

COMPOSITION

Reagent: R

Good's buffer, pH 6.15

CNP-G₃ 2.27 mmol/L

Sodium azide < 0.1 % (w/w)

Also contains calcium and chloride salts, as well as potassium thiocyanate for optimal performance.

MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2

- CONT-0060 ELITROL I

- CONT-0160 ELITROL II

- Normal saline solution (NaCl 9 g/L).

- Analyzers.

- General Laboratory equipment (e.g. pipette).

- Do not use materials that are not required as indicated above.

PRECAUTIONS FOR USE AND WARNINGS

- Consult Safety Data Sheet (SDS) for a proper handling.

- Reagent R contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.

- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contamination.

- Saliva and sweat contain amylase. It is therefore recommended to wear gloves and a mask to avoid the contamination of the reagent.

- Contact with acids liberates very toxic gas.

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