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**English - EN**

**INTENDED USE**

*For Total Bilirubin* : BILIRUBIN TOTAL 4+1 is an *in vitro* diagnostic reagent intended for the quantitative determination of total bilirubin in human serum and plasma samples on adults and children over 10 days on analyzers or semi-automatic analyzers.

*For Direct Bilirubin* : BILIRUBIN DIRECT 4+1/ BILIRUBIN DIRECT are *in vitro* diagnostic reagents intended for the quantitative determination of direct bilirubin in human serum and plasma samples on analyzers or semi-automatic analyzers.

These *in vitro* diagnostic devices are for professional use only.

**CLINICAL SIGNIFICANCE** <sup>(1-2)</sup>

Approximately 80-85 % of the bilirubin originates from the degradation of hemoglobin's heme moiety, is subsequently transported into the liver where it is rapidly rendered soluble by glucuronidation, is excreted into biliary canaliculi, and is finally hydrolyzed in the gastrointestinal tract.

Unconjugated bilirubin serum concentration increases in case of overproduction of bilirubin (acute or chronic haemolytic anemias) and in case of disorders of bilirubin metabolism or transport (impaired uptake by liver cells: Gilbert's syndrome; defects in the conjugation reaction: Crigler-Najjar syndrome). Reduced excretion (hepatocellular damage such as hepatitis or cirrhosis; Dubin-Johnson and Rotor syndromes) and obstruction to the flow of bile (most often produced by gallstones or by tumors) induce an important elevation of conjugated bilirubin and in a minor extent an increase of unconjugated bilirubin (conjugated hyperbilirubinemia). In clinical practice, serum bilirubin measurement in children and adults is indicated to help screen and follow-up of liver disorders and hemolytic anemia.

**LIMITATION OF USE**

The quantitative assay of total & direct bilirubin alone cannot be used to diagnose a disease or a specific pathology.

The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

**METHOD & PRINCIPLE** <sup>(2)</sup>

Malloy-Evelyn modified. End point.

In the presence of accelerator (cetrimide, Total Bilirubin), conjugated and unconjugated bilirubin react with diazotized sulphanic acid to form azobilirubin. In the absence of accelerator (Direct Bilirubin), only conjugated bilirubin reacts.

Sulphanilic acid + NaNO <sub>2</sub>	→	Diazotized sulphanic acid
Bilirubin + Diazotized sulphanic acid	→	Azobilirubin

The increase of absorbance at 546 nm is proportional to bilirubin concentration.

**COMPOSITION**

**Bilirubin Total 4+1 :**

<b>Reagent 1: R1</b>	29 mmol/L
Sulphanilic acid	29 mmol/L
Cetrimide	29 mmol/L

**Bilirubin Direct 4+1/Bilirubin Direct :**

<b>Reagent 1: R1</b>	29 mmol/L
Sulphanilic acid	29 mmol/L

**Bilirubin Total & Direct 4+1/Bilirubin Direct :**

<b>Reagent 2: R2</b>	11 mmol/L
Sodium nitrite	11 mmol/L

**MATERIALS REQUIRED BUT NOT PROVIDED**

- CALI-0550 ELICAL 2

- CONT-0060 ELITROL I

- CONT-0160 ELITROL II

- Normal saline solution (NaCl 9 g/L).

- Biochemistry analyzer equipped with required filters (Refer to § PROCEDURE) or semi-automatic analyzers.

- General Laboratory equipment (e.g. pipette).

- Do not use materials that are not required as indicated above.

**PRECAUTIONS FOR USE AND WARNINGS**

- The reagents R1 contain sulphanilic acid . They may produce an allergic reaction.

- Consult Safety Data Sheet (SDS) for a proper handling.

- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contamination.

- Do not interchange reagent vials from different kits.

**STABILITY**

**Store at 2-8 °C and protect from light. Do not freeze.** Do not use after expiration dates indicated on the vial labels.

*On board stability :*

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

**PREPARATION**

The reagents are ready to use.

**PRODUCT DETERIORATION**

- Reagent R1 of Bilirubin Total 4+1 can be slightly cloudy with no impact on product performances.

- The reagent R1 of Bilirubin Direct 4+1/Bilirubin Direct and the reagent R2 of Bilirubin Total & Direct 4+1/ Bilirubin Direct should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).

- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

**SAMPLES**

**Specimen** <sup>(3)</sup>

- Serum.

- Plasma (lithium heparin).

- Using any other specimen type should be validated by the laboratory.

**Warnings and precautions**

- Protect the samples from light before and during analysis<sup>(9)</sup>

- Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

**Storage and stability** <sup>(3)</sup>

- 1 day (Total bilirubin) and 2 days (Direct bilirubin) at room temperature.

- 7 days at 2-8°C.

- 6 months at -20°C.

**REFERENCE VALUES** <sup>(4,5)</sup>

<i>Serum/plasma</i>	mg/dL	µmol/L
<i>Total bilirubin</i>	≤ 1.2	≤ 21.0
Adults and children over 10 days old		

Total bilirubin concentrations observed for breastfed children remain more elevated for a longer period.

*Direct bilirubin*

Adults	≤ 0.2	≤ 3.4
--------	-------	-------

*Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

**INSTALLATION AND USE**

*For use on Selectra Pro analyzers :*

- Consult operator manual.

- **Special Programming instructions: Programming special instructions is mandatory when some combinations of tests are performed together on the analyzer.** Refer to Instructions For Use of ACID SOLUTION & SYSTEM CLEANING SOLUTION for adequate programming (See PIT-SOL).

**PROCEDURE**

**A) Total Bilirubin**

*Manual Procedure*

*Wavelength : 550 nm (λ. main) - 700 nm (λ. sub)*

*Optical path : 1 cm*

*Sample/reagent ratio : 1:20*

*Temperature : 37 °C*

*Read against distilled water.*

	CALIBRATION	TEST
<b>Reagent R1</b>	800 µL	800 µL
<b>Callibrator</b>	50 µL	-
<b>Sample</b>	-	50 µL

Mix and read the absorbances A<sub>1</sub> after an incubation of 5 minutes

$\Delta A_1 = A_{\lambda_1(550\text{ nm})} - A_{\lambda_1(700\text{ nm})}$

<b>Reagent R2</b>	200 µL	200 µL
Mix and read the absorbances A <sub>2</sub> after an incubation of 5 minutes		
	$\Delta A_2 = A_{\lambda_2(550\text{ nm})} - A_{\lambda_2(700\text{ nm})}$	

*Automatic Procedure*

These reagents may be used on several automatic analyzers. For Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

For Selectra ProXS users, a 700 nm additional filter is required.

**B) Direct bilirubin**

*Manual Procedure*

*Wavelength : 550 nm*

*Optical path : 1 cm*

*Sample/reagent ratio : 1:10*

*Temperature : 37 °C*

*Read against distilled water.*

	CALIBRATION	TEST
<b>Reagent R1</b>	800 µL	800 µL
<b>Calibrator</b>	100 µL	-
<b>Sample</b>	-	100 µL
Mix and read the absorbance (A <sub>1</sub> ) after an incubation of 5 minutes		
<b>Reagent R2</b>	200 µL	200 µL
Mix and read the absorbance (A <sub>2</sub> ) after an incubation of 50 seconds		

*Automatic Procedure*

These reagents may be used on several automatic analyzers. For Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

**In the Selectra Pro application, the offset must be set to :**

**- 0.05 mg/dL (- 0.9 µmol/L).**

**CALCULATION**

**A) Total Bilirubin**

$[\Delta A_2 - \Delta A_1 \times Fdil]$	<i>Sample</i>	x n
$[\Delta A_2 - \Delta A_1 \times Fdil]$	<i>Calibrator</i>	

Fdil = Dilution factor = 0.81 (manual procedure)

n = Calibrator concentration

**B) Direct Bilirubin**

$[A_2 - A_1 \times Fdil]$	<i>Sample</i>	x n
$[A_2 - A_1 \times Fdil]$	<i>Calibrator</i>	

Fdil = Dilution factor = 0.82 (manual procedure)

n = Calibrator concentration

**Conversion factor:** mg/dL x 17.1 = µmol/L

**CALIBRATION**

ELICAL 2 is traceable to SRM 916a reference material.

**Calibration frequency :** The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

**QUALITY CONTROL**

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay.

Controls have to be performed :

- prior to assaying patient samples,

- at least once per day,

- after every calibration,

- and/or in accordance with laboratory and regulatory requirements.

Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

**WASTE MANAGEMENT**

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

**PERFORMANCES**

Performances were obtained on Selectra ProM, following CLSI technical recommendations, under controlled environmental conditions.

- **Measuring range**

a) *Total bilirubin*

0.25 - 25.00 mg/dL (4.3 - 427.6 µmol/L)

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 60.00 mg/dL (1026.3 µmol/L).

Do not report results outside this extended range.

b) *Direct bilirubin*

0.08 - 10.55 mg/dL (1.4 - 180.4 µmol/L).

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 50.00 mg/dL (855.2 µmol/L).

Do not report results outside this extended range.

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

**Limit of Detection (LoD) and Limit of Quantification (LoQ)**

a) *Total bilirubin*

LoD = 0.04 mg/dL (0.7 µmol/L)

LoQ = 0.15 mg/dL (2.6 µmol/L)

b) *Direct bilirubin*

LoD = 0.01 mg/dL (0.2 µmol/L)

LoQ = 0.08 mg/dL (1.4 µmol/L)

- **Precision**

Imprecision data has been obtained on 2 Selectra ProM analyzers over 20 days (2 runs per day, tests performed in duplicate). Representative results are presented below.

a) *Total bilirubin*

		Mean	Within-run	Total	
	n	mg/dL	µmol/L	CV (%)	
<b>Level 1</b>	80	1.15	19.7	1.8	5.0
<b>Level 2</b>	80	4.08	69.8	0.4	3.1
<b>Level 3</b>	80	14.61	249.9	0.5	2.9

b) *Direct bilirubin*

		Mean	Within-run	Total	
	n	mg/dL	µmol/L	CV (%)	
<b>Level 1</b>	80	0.36	6.2	3.8	5.2
<b>Level 2</b>	80	1.51	25.8	1.9	5.3
<b>Level 3</b>	80	3.99	68.2	0.9	4.7

- **Correlation**

a) *Total bilirubin*

A comparative study has been performed between BILIRUBIN TOTAL 4+1 reagent on a Selectra ProM analyzer and a similar commercially available system on 100 human serum samples.

The sample concentrations ranged from 0.32 to 23.02 mg/dL (5.5 - 393.7 µmol/L).

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## BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFÍA/BIBLIOGRAFIA

- Higgins, T., et al., *Hemoglobin, Iron, and Bilirubin, Tietz Fundamentals of Clinical Chemistry*, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R. Bruns, D.E. (W.B. Saunders eds.), (2008), 509.
- Dufour, D.R. *The liver: Function and Chemical Pathology, Clinical Chemistry: Theory, Analysis, Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds.), (2010), 586 and appendix.
- Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
- Wu, A.H.B., *Clinical guide to laboratory tests*, 4<sup>th</sup> Ed., (W.B. Saunders eds.), (2006), 172.
- Howard, C.R., *Breastmilk & Breastfeeding Jaundice, Pediatric Clinical Advisor: Instant Diagnosis and Treatment*, 2<sup>nd</sup> Ed., Garfunkel, L.C., Kaczorowski, J.M. & Christy, C., (Mosby Inc. eds.), (2007), 82.
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), 59, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press (1997).
- Young D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press (1995).

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- Symbols used on our documentation are defined on ISO-15223-1 standard, except for some presented in the symbols glossary available on the VitalScientific Website. (Symbols glossary).

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- Special Programming instructions : see § INSTALLATION AND USE

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