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BIBLIOGRAPHIE/BIBLIOGRAPHY
BIBLIOGRAFÍA/BIBLIOGRAFIA

- Rifai, N., Warnick, G.R., Remaley, A.T., *Lipids, lipoproteins, apolipoproteins and other cardiovascular risk factors*. Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.),(2008), 402.
- Burnett, J.R., *Coronary Artery Disease: Lipid metabolism*. Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds.), (2010), 691 and appendix.
- Langlois, M.R., *et al.* for the European Atherosclerosis Society (EAS) and the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) Joint Consensus Initiative, Clin. Chem. Lab. Med., (2020), **58**, 496.
- Allain, C.C., et al., Clin. Chem., (1974), 20, 470.
- Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/ LAB/99.1 Rev.2.
- Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACCPress, (1997).
- Young D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACCPress, (1995).

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English - EN

INTENDED USE

CHOLESTEROL SL is an *in vitro* diagnostic reagent intended for the quantitative determination of total cholesterol in human serum and plasma samples on analyzers or semi-automatic analyzers. The standard is intended for the calibration of reagent. These *in vitro* diagnostic devices are for professional use only.

CLINICAL SIGNIFICANCE ⁽¹⁻³⁾

Total cholesterol in the serum is derived from dietary sources or is synthesized endogenously, mainly in hepatic and intestinal cells. Cholesterol is an important structural component of cell and organelle membranes. It is a precursor of bile acids, vitamin D and steroid hormones. Cholesterol, being insoluble in water, circulates in association with lipoproteins (HDL, LDL, VLDL and chylomicrons). In practice, total cholesterol measurement is requested to assess predisposition of patients to cardiovascular risk as part of a lipid profile and to monitor associated therapeutic strategies. Total Cholesterol measurement is also important to help diagnose hyperlipoproteinemia.

LIMITATION OF USE

The quantitative assay of total cholesterol alone cannot be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.