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#### BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFIA/BIBLIOGRAFIA

- Kaplan, J.M., First, M.R., *Renal Function. Clinical Chemistry: Theory, Analysis, Correlation*. 5<sup>th</sup> Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 567 and appendix.
- Delaney, M.P., Price, C.P., Lamb, E.J., *Kidney Function and Disease. Tietz Fundamentals of Clinical Chemistry*. 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 631.
- Wu, A. H. B., *Clinical guide to laboratory tests*, 4<sup>th</sup> Ed., (W.B. Saunders eds.), (2006), 316.
- Vasiliades, J., Clin. Chem., (1976), **22**, 1664.
- Guder, W.G., *et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/ LAB/99.1 Rev.2.
- Ceriotti, F. *et al., Reference Intervals for Serum Creatinine Concentrations: Assessment of Available Data for Global Application, Clin\_Chem.*, (2008), **54**, 559.
- Young D.S., **Effects of drugs on clinical laboratory tests**, 4<sup>th</sup> Ed., AACCC Press (1995).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), 59, 263.
- Young, D.S., **Effects of preanalytical variables on clinical laboratory tests**, 2<sup>nd</sup> Ed., AACCC Press (1997).

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- Symbols used on our documentation are defined on ISO-15223-1 standard, except for some presented in the symbols glossary available on the VitalScientific Website. (Symbols glossary).

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### English - EN

#### INTENDED USE

CREATININE JAFFE is an *in vitro* diagnostic reagent intended for the quantitative determination of creatinine in human serum samples on analyzers or semi-automatic analyzers.

The standard is intended for the calibration of the reagent.

These *in vitro* diagnostic devices are for professional use only.

#### CLINICAL SIGNIFICANCE <sup>(1-3)</sup>

Creatinine is the waste product of creatine, a compound present mainly in muscle tissues. It is usually produced a relatively constant rate and is freely filtered by the glomerulus, which makes it an excellent marker of renal function, with serum creatinine concentration increasing in case of acute or chronic renal insufficiency.

In practice, creatinine measurement in blood performed to help diagnose kidney diseases, to screen for early kidney failure for patients at risk, to monitor treatments against these pathologies or to monitor renal function during treatments that may impair this function.

#### LIMITATION OF USE

The quantitative assay of creatinine alone can not be used to diagnose a disease or a specific pathology.

The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

#### METHOD & PRINCIPLE <sup>(4)</sup>

Jaffe – Kinetic.

In alkaline medium, creatinine reacts with picrate to form a red complex. The rate of increase of absorbance measured at 505 nm is proportional to creatinine concentration.

#### COMPOSITION

**Reagent 1: R1**  
Sodium hydroxide 0.31 mol/L  
Also contains phosphate

**Reagent 2: R2**  
Picric acid 8.7 mmol/L  
**Standard: Std**  
Creatinine 2 mg/dL  
177 µmol/L

#### MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2  
- CONT-0060 ELITROL I  
- CONT-0160 ELITROL II  
- Normal saline solution (NaCl 9 g/L).  
- Analyzers or semi-automatic analyzers.  
- General Laboratory equipment (e.g. pipette).  
- Do not use materials that are not required as indicated above.

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## PRECAUTIONS FOR USE AND WARNINGS

- The reagent R1 is classified as hazardous :

- Consult Safety Data Sheet (SDS) for a proper handling.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Do not interchange reagent/standard vials from different kits.

### STABILITY

**Store at 2-25 °C and protect from light. Do not freeze.**

Do not use after expiration dates indicated on the vial labels.

The standard should be immediately and tightly capped to prevent contamination and evaporation.

**On board stability :**

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

### PREPARATION

The reagent and standard are ready to use.

### PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).

- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

### SAMPLES

#### Specimen <sup>(5)</sup>

- Serum.

- Using any other specimen type should be validated by the laboratory.

#### Warnings and precautions

Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

#### Storage and stability <sup>(5)</sup>

- 7 days at 20-25 °C.

- 7 days at 4-8 °C.

- 3 months at - 20 °C.

### REFERENCE VALUES <sup>(3,6)</sup>

<i>Serum</i>	mg/dL	µmol/L
Adults		
Men <span> </span> :	0.8 - 1.3	71 - 115
Women <span> </span> :	0.6 - 1.2	53 - 106

***Note :*** *The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

#### ❖INSTALLATION AND USE

**For use on Selectra Pro analyzers :**

- Consult operator manual.

- **Special Programming Instructions: Programming special instructions is mandatory when some combinations of tests are performed together on the analyzer.** Refer to Instructions For Use of ACID SOLUTION & SYSTEM CLEANING SOLUTION for adequate programming (See PIT-SOL).

### PROCEDURE

*Manual Procedure*

<i>Wavelength<span> </span>:</i>	<i>505 nm</i>
<i>Optical path<span> </span>:</i>	<i>1 cm</i>
<i>Sample/reagent ratio<span> </span>:</i>	<i>1:10</i>
<i>Temperature<span> </span>:</i>	<i>37<span> </span>°C</i>
<i>Read against distilled water.</i>	

<b>Working reagent (1 volume of R1 + 1 volume of R2)</b>	1000 µL
<b>Sample or Standard</b>	100 µL

*Mix and read the absorbance (A1) 10 seconds after the sample or standard addition. Exactly 2 minutes after the first reading, take second reading (A2).*

#### Automatic Procedure

These reagents may be used in several automatic analyzers. For Selectra Analyzers, validated applications are available on request.

### CALCULATION

(A2 - A1) Sample x n = standard concentration
(A2 - A1) Standard

**Conversion factor :** mg/dL x 88.40 = µmol/L
mg/dL x 10 = mg/L

### CALIBRATION

Creatinine Standard 2 mg/dL and ELICAL 2 are traceable to ID-MS (Isotope Dilution - Mass Spectrometry) reference method.

**Calibration frequency :** The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

### QUALITY CONTROL

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay.

Controls have to be performed :

- prior to assaying patient samples,

- at least once per day,

- after every calibration,

- and/or in accordance with laboratory and regulatory requirements.

Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

## WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

## PERFORMANCES

Performances were obtained on Selectra Pro, following CLSI technical recommendations, under controlled environmental conditions.

### - Measuring range

0.20 - 15.00 mg/dL (18 - 1326 µmol/L)

Do not report results outside the measuring range.

**Limit of Detection (LoD) and Limit of Quantification (LoQ)**

**LoD** = 0.05 mg/dL (4 µmol/L)

**LoQ** = 0.20 mg/dL (18 µmol/L)

### - Precision

Imprecision data has been obtained on 2 Selectra Pro analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented below.

		Mean	Within-run	Total	
	n	mg/dL	µmol/L	CV (%)	
<b>Level 1</b>	80	0.95	84	1.1	1.8
<b>Level 2</b>	80	1.51	133	1.1	1.7
<b>Level 3</b>	80	4.98	440	0.9	1.6

### - Correlation

A comparative study has been performed between CREATININE JAFFE reagent on a Selectra ProM analyzer and a similar commercially available system on 108 human serum samples.

The sample concentrations ranged from 0.16 to 14.80 mg/dL (14 - 1308 µmol/L).

The results are as follows :

Correlation coefficient : (r) = 0.999

Linear regression: y = 0.939x + 0.14 mg/dL (12 µmol/L).

### - Limitations/Interferences

- Antibiotics containing cephalosporin lead to significant false-positive values.<sup>(7)</sup>

- The presence of ketone bodies can cause artificially high results.<sup>(1)</sup>

- Studies have been performed to determine the level of interference from different compounds. The following creatinine levels were tested: 1.50 mg/dL and 5.00 mg/dL.

No significant interference is defined by a recovery ±10% of the initial value.

**Triglycerides:** No significant interference up to 2 100 mg/dL (23.7 mmol/L).

**Unconjugated bilirubin:** No significant interference up to 15.0 mg/dL (257 µmol/L).

**Conjugated bilirubin:** No significant interference up to 5.9 mg/dL (101 µmol/L).

**Hemoglobin:** No significant interference up to 500 mg/dL.

**Glucose :** No significant interference up to 500 mg/dL (27.8 mmol/L).

**Creatine:** No significant interference up to 10.0 mg/dL (762 µmol/L).

**Uric acid :** No significant interference up to 20.0 mg/dL (1190 µmol/L).

**L-dopa :** No significant interference up to 20.0 mg/dL.

**Áscorbic Ácido:** No significant interference up to 20 mg/dL.

**Acetaminophen :** No significant interference up to 30 mg/dL.

**Acetylsalicylic acid :** No significant interference up to 200 mg/dL.

**Methyl dopa :** No significant interference up to 2 mg/dL.

**Pyruvate :** No significant interference up to 3.0 mg/dL (341 µmol/L)

**Proteins :** No significant interference between 3.0 g/dL and 11.0 g/dL.

- Do not use icteric samples.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.<sup>(8)</sup>

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.<sup>(7,9)</sup>

#### - On board stability/Calibration frequency

**On Board Stability:** 2 days

**Calibration frequency:** At the end of onboard stability period,change the bottle with a new one filled with fresh reagent.

Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

*These performances have been obtained using Selectra Pro analyzer. Results may vary if a different instrument or a manual procedure is used.*

*❖The performances of applications not validated by VitalScientific are not warranted and must be defined by the user.*

## DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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