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2. Burnett, J.R., *Coronary Artery Disease: Lipid Metabolism. Clinical Chemistry: Theory, Analysis, Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A. Pesce, A.J., (Mosby Inc. eds.), (2010), 691 and appendix.

3. Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), *Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP)*, *JAMA*, (2001), **285**, 2486.

4. Fossati, P. & Prencipe, L., *Clin. Chem.*, (1982), **28**, 2077.

5. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.

6. Langlois, M.R., *et al.* for the European Atherosclerosis Society (EAS) and the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) Joint Consensus Initiative, *Clin. Chem. Lab. Med.*, (2020), **58**, 496.

7. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.

8. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACCC Press (1997).

9. Young D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACCC Press (1995).

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- Symbols used on our documentation are defined on ISO-15223-1 standard, except for some presented in the symbols glossary available on the VitalScientific Website. (Symbols glossary).

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
- **Special Programming instructions :** see § INSTALLATION AND USE

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## English - EN

### INTENDED USE

TRIGLYCERIDES SL and TRIGLYCERIDES MONO SL NEW are an *in vitro* diagnostic reagents intended for the quantitative determination of triglycerides in human serum and plasma samples on analyzers or semi-automatic analyzers.

The standard is intended for the calibration of the reagent.

These *in vitro* diagnostic devices are for professional use only.

### CLINICAL SIGNIFICANCE <sup>(1-3)</sup>

Triglycerides constitute 95% of tissue storage fat and their main role is to provide energy for the cell. They are synthesized both in the intestine from dietary fats and in liver from dietary carbohydrates, and are then transported in blood by chylomicrons and VLDL. Situations where high levels of triglycerides are observed can be obesity and overweight, physical inactivity, cigarette smoking, excess alcohol intake, high carbohydrate diets, several diseases such as type 2 diabetes, chronic renal failure, nephrotic syndrome, certain drugs and genetic dyslipoproteinemia.

In practice, triglyceride measurement is requested to assess predisposition of patients to cardiovascular risk as part of a lipid profile and to monitor associated therapeutic strategies.

### LIMITATION OF USE

The quantitative assay of triglycerides alone cannot be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

### METHOD & PRINCIPLE <sup>(4)</sup>

Enzymatic / PAP - End point.

*Lipoprotein Lipase* → Triglycerides + H<sub>2</sub>O → Glycerol + Fatty acids

*Glycerol kinase* → Glycerol + ATP → Glycerol-3-Phosphate+ ADP

*GPO* → Glycerol-3-Phosphate + O<sub>2</sub> → DHA-P + H<sub>2</sub>O<sub>2</sub>

*Peroxidase* → H<sub>2</sub>O<sub>2</sub> + 4-AAP + p-Chlorophenol → Quinoneimina

GPO = Glycerol-3-phosphate oxidase
DHA-P = Dihydroxyacetone-Phosphate
4-AAP = Amino-4-antipirina

### COMPOSITION

**Reagent : R**

Good's buffer, pH 7.00

p-Chlorophenol	2.7	mmol/L
ATP	3.15	mmol/L
Amino-4-antipyrine (4-AAP)	0.31	mmol/L
Lipoprotein lipase	≥ 2 000	U/L
Glycerol kinase	≥ 500	U/L
Glycerol-3-phosphate oxidase	≥ 4 000	U/L
Peroxidase	≥ 500	U/L
Sodium azide	< 0.1	% (w/w)
Also contains magnesium salts, FAD and surfactants for optimal performance		
<b>Standard: Std</b> (Ref <span> </span> : TGM.L-0427/497/0517/0707)		
Glycerol (triglycerides equivalent)	200	mg/dL
	2.26	mmol/L
Sodium azide	<	0.1 <span> </span> % (w/w)

### MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2

- CONT-0060 ELITROL I

- CONT-0160 ELITROL II

- Normal saline solution (NaCl 9 g/L).

- Analyzers or semi-automatic analyzers.

- General Laboratory equipment (e.g. pipette).

- Do not use materials that are not required as indicated above.

### PRECAUTIONS FOR USE AND WARNINGS

- Consult Safety Data Sheet (SDS) for a proper handling.

- Reagent R and standard Std contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.

- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contamination.

### STABILITY

**Store at 2-8 °C and protect from light. Do not freeze.** Do not use after expiration dates indicated on the vial labels.

The standard should be immediately and tightly capped to prevent contamination and evaporation.

**On board stability :**

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

### PREPARATION

The reagent and standard are ready to use.

### PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).

- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

### SAMPLES

**Specimen <sup>(2,5)</sup>**

- Serum.

- Plasma (lithium heparin).

- Using any other specimen type should be validated by the laboratory.

**Warnings and precautions**

- A non-fasting or fasting sample is suitable for the determination of a lipid profile. A repeat of lipid profile on a fasting sample could be performed in cases where a nonfasting Triglyceride result> 400 mg/dL (4.5 mmol/L) or hypertriglyceridemia is known <sup>(6)</sup>

- Collect the sample in tubes and stoppers free of glycerol.<sup>(2)</sup>

- Separate from cells within 2 hours.<sup>(2)</sup>

- Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

**Storage and stability <sup>(6)</sup>**

- 2 days at room temperature

- 7 days at 2-8°C

- 1 year at -20°C

### REFERENCE VALUES <sup>(6)</sup>

Most recent publications recommend adapting triglyceride thresholds as part of an overall risk assessment. At the laboratory level, the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) recommends that the following concentrations be reported as abnormal:

<i>Serum/plasma</i>	mg/dL	mmol/L
Fasting sample	≥ 150	≥ 1.7
Non-fasting sample	≥ 175	≥ 2.0

***Note:** Laboratories should follow locally applicable recommendations for lipid targets if they differ from those reported above.*

### INSTALLATION AND USE

**For use on Selectra Pro analyzers :**

- Consult operator manual.

- **Special Programming instructions: Programming special instructions is mandatory when some combinations of tests are performed together on the analyzer.** Refer to Instructions For Use of ACID SOLUTION & SYSTEM CLEANING SOLUTION for adequate programming (See PIT-SOL).

☛In spite of these instructions, TRIGLYCERIDES SL / TRIGLYCERIDES MONO SL NEW reagent may be prone to on-board contamination; repeat testing when results are not as expected.

### PROCEDURE

*Manual Procedure*

*Wavelength :* 505 nm

*Optical path :* 1 cm

*Sample/ Reagent ratio :* 1:100

*Temperature:* 37 °C

*Read against reagent blank.*

	CALIBRATION	TEST
<b>Reagent R</b>	1 000 µL	1 000 µL
<b>Standard/ Calibrator</b>	10 µL	-
<b>Sample</b>	-	10 µL

*Mix and read the absorbances (A) after an incubation of 10 minutes.*

### Automatic Procedure

These reagents may be used on several automatic analyzers. For Selectra Analyzers, validated applications are available on request. (Si code barre dans l'IFU) For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

### CALCULATION

A *Sample* × *n* = *n* = Calibrator/ standard
A *Calibrator/ Standard* × *concentration*

**Conversion factor :** mg/dL x 0.0113 = mmol/L
mg/dL x 0.01 = g/L

### CALIBRATION

ELICAL 2 and Triglycerides Standard 200 mg/dL are traceable to ID-GC-MS (Isotope Dilution - Gas Chromatography - Mass Spectrometry) reference method.

**Calibration frequency :** The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

### QUALITY CONTROL

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay.

Controls have to be performed :
- prior to assaying patient samples,
- at least once per day,
- after every calibration,
- and/or in accordance with laboratory and regulatory requirements.

Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

### WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

### PERFORMANCES

Performances were obtained on Selectra ProM, following CLSI technical recommendations, under controlled environmental conditions.

### - Measuring range

30 - 1 000 mg/dL (0.34 - 11.30 mmol/L)

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 5 000 mg/dL (56.50 mmol/L).

Do not report results outside this extended range.

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

**Limit of Detection (LoD) and Limit of Quantification (LoQ)**

LoD = 3 mg/dL (0.03 mmol/L)

LoQ = 10 mg/dL (0.11 mmol/L)

### - Precision

Imprecision data has been obtained on 2 Selectra ProM analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented below.

		Mean	Within-run	Total	
	n	mg/dL	mmol/L	CV (%)	
<b>Level 1</b>	80	44	0.50	2.0	3.8
<b>Level 2</b>	80	131	1.48	0.9	2.3
<b>Level 3</b>	80	267	3.02	1.2	2.4

### - Correlation

A comparative study has been performed between TRIGLYCERIDES SL reagent on a Selectra ProM analyzer and a similar commercially available system on 99 human serum samples. The sample concentrations ranged from 30 to 957 mg/dL (0.34 - 10.81 mmol/L).

The results are as follows :

Correlation coefficient : (r) = 0.999

Linear regression: y = 1.019 x + 1 mg/dL (0.01 mmol/L)

### - Limitations/Interferences

Studies have been performed to determine the level of interference from different compounds.

The following triglycerides levels were tested: 133 mg/dL and 266 mg/dL.

No significant interference is defined by a recovery ≤±10% of the initial value.

**Unconjugated Bilirubin**: No significant interference up to 15.0 mg/dL (257 µmol/L).

**Conjugated Bilirubin**: No significant interference up to 5.9 mg/dL (101 µmol/L).

**Hemoglobin** : No significant interference up to 125 mg/dL.

**Uric acid** : No significant interference up to 24.2 mg/dL (1440 µmol/L).

**Ascorbic acid** : No significant interference up to 2.0 mg/dL. Concentrations above the therapeutic levels will interfere and cause erroneous results.

**Methyl-dopa** : No significant interference up to 1.0 mg/dL.

Do not use hemolyzed or icteric samples.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.<sup>(7)</sup>

- Results can be falsely lowered by significant levels in the sample of NAC (*N*-Acetyl-Cysteine), NAPQI (metabolite of acetaminophene (paracetamol)) or metazolole.

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.<sup>(8-9)</sup>

#### - On board stability/Calibration frequency

**On Board Stability:** 28 days

**Calibration frequency:** 14 days

Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

*These performances have been obtained using Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used.*

☛*The performances of applications not validated by VitalScientific are not warranted and must be defined by the user.*

### DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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