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

Enzimática / PAP - Ponto final.

EHSPT = N-etil-N-(2-hidroxi-3-sulfopropil)-m-toluidina
4-AAP = amino-4-antipirina

COMPOSIÇÃO

Reagente : R

Tampão pH 7.00 (20-25°C)

EHSPT	0.72	mmol/L
Amino-4-antipirina	0.37	mmol/L
Uricase	≥ 150	U/L
Peroxidase	≥ 12 000	U/L
Azida de sódio	< 0.1	% (p/p)

Padrão: Std (Ref : AUML-0427/0497/0507/0707)

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English - EN

INTENDED USE

URIC ACID MONO SL is an *in vitro* diagnostic reagent intended for the quantitative determination of Uric Acid in human serum, plasma and urine samples on analyzers or semi-automatic analyzers.

The standard is intended for the calibration of reagent. These *in vitro* diagnostic devices are for professional use only.

CLINICAL SIGNIFICANCE ⁽¹⁻³⁾

Uric acid is the major product of the catabolism of endogenous and exogenous purines (adenosine and guanosine). Uric acid is poorly soluble in water. Therefore, in case of high serum concentration, urate crystals can form and deposit in joints, triggering painful inflammations (gout), or can damage kidneys. Increased serum uric acid level can be caused either by increased production (increased intake of purines, increased nucleic acid turn-over particularly in case of certain cancers or following anti-cancer treatments, genetic metabolic disorders such as Lesch-Nyhan syndrome, psoriasis) or by decreased excretion (renal failure, drugs such as diuretics). In the case of pre-eclampsia, serum uric acid can be increased due to both mechanisms. Decreased serum uric acid is more uncommon. It can occur for example in impaired renal elimination such as in Fanconi syndrome or in Hodgkin's disease.

When uric acid concentration is abnormally high in urines, there is a risk of stone formation.

☛In practice, the determination of uric acid in serum is indicated to help diagnose pathologies involving hyperuricemia such as gout. The determination of uric acid in urines is indicated to identify the nature of kidney stones and prevent those from recurring.

LIMITATION OF USE

The quantitative assay of Uric Acid alone cannot be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

METHOD & PRINCIPLE ⁽⁴⁾

Enzymatic / PAP - End Point

Uricase

{\displaystyle \rightarrow }

Uric acid + 2H₂O + O₂ → Allantoïne + CO₂ + H₂O₂

Peroxidase

{\displaystyle \rightarrow }

2H₂O₂ + 4-AAP + EHSPT → Quinoneimine + 4H₂O

EHSPT = N-Ethyl-N-(2-Hydroxy-3-Sulfopropyl) *m*-Toluidine
4-AAP = Amino-4-antipyrine

COMPOSITION

Reagent: R

Buffer pH 7.00 (20-25°C)

EHSPT	0.72	mmol/L
Amino-4-antipyrine	0.37	mmol/L
Uricase	≥ 150	U/L
Peroxidase	≥ 12 000	U/L
Sodium azide	< 0.1	% (w/w)
Standard: Std (Ref : AUML-0427/0497/0507/0707)		
Uric acid	6	mg/dL
	357	µmol/L
Sodium azide	<	0.5 % (w/w)

MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- Normal saline solution (NaCl 9 g/L).
- Analyzers or semi-automatic analyzers.
- General Laboratory equipment (e.g. pipette).
- Do not use materials that are not required as indicated above.

PRECAUTIONS FOR USE AND WARNINGS

- The reagent R is classified as hazardous :

DANGER. May damage fertility. May damage the unborn child. Wear protective gloves/protective clothing/eye protection/face protection. Do not handle until all safety precautions have been read and understood. IF exposed or concerned: Get medical advice/attention.

Obtain Safety data sheet (SDS) before use for a proper handling.

- Reagent R and Standard Std contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.

- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.

STABILITY
Store at 2-8 °C and protect from light. Do not freeze. Do not use after expiration dates indicated on the vial labels.

The standard should be immediately and tightly capped to prevent contamination and evaporation.

On board stability :

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagent and standard are ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).

- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

SAMPLES Specimen ⁽²⁾

- Serum

- Plasma (lithium heparin)

- Urine

- Using any other specimen type should be validated by the laboratory.

Warnings and precautions

- To prevent urate precipitation, urine samples may be adjusted to pH>8.0 with NaOH.⁽²⁾

- Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

Storage and stability ⁽²⁾

Serum/ plasma

- 3-5 days at 2-8°C

- 6 months at -20°C

Urine (Alkalinized)

- 3 days at room temperature

Do not refrigerate urine samples

REFERENCE VALUES ⁽¹⁾

<i>Serum/plasma</i>	mg/dL	µmol/L		
Men	3.5 - 7.2	208 - 428		
Women	2.6 - 6.0	155 - 357		
<i>Urine (24 h collection)</i>	mg/24h	mmol/24		
	250 - 750	1.48 - 4.43		
			For a urinary volume of 1.5 L/24h	
	mg/dL	mmol/L		
	16.7 - 50	0.99 - 2.97		

With a purine-free diet, excretion may decrease by 20 to 25%

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

INSTALLATION AND USE

For use on Selectra Pro analyzers :

- Consult operator manual.

- **Special Programming instructions: Programming special instructions is mandatory when some combinations of tests are performed together on the analyzer.** Refer to Instructions For Use of ACID SOLUTION & SYSTEM CLEANING SOLUTION for adequate programming (See PIT-SOL).

PROCEDURE

Manual Procedure

Wavelength : 546 nm
Optical path : 1 cm
Sampler/ Reagent ratio : 1:40
Temperature : 37 °C
Urine samples must be diluted 1:10 with NaCl 9 g/L solution before measurement.
Read against reagent blank.

	CALIBRATION	TEST
Reagent R	1 000 µL	1 000 µL
Standard/ Calibrator	25 µL	-
Sample	-	25 µL

Mix and read the absorbances (A) after an incubation of 5 minutes.

Automatic Procedure

These reagents may be used on several automatic analyzers. For Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

Urine samples must be diluted 1:10 with NaCl 9 g/L solution before measurement. For users of Selectra TouchPro software, urine dilution is performed automatically.

CALCULATION

A
Sample

x
n
=
Calibrator/
standard
concentration

{\displaystyle A_{Sample}\times n={\frac {Calibrator}{standard\ concentration }}}

For the calculation of uric acid concentration in urine, multiply the result by the dilution factor (10). For users of Selectra TouchPro software, the results take the dilution factor into account.

Conversion factor : mg/dL x 59.48 = µmol/L
mg/dL x 0.059 = mmol/L

CALIBRATION

ELICAL 2 and Uric Acid Standard 6 mg/dL are traceable to ID-MS (Isotope Dilution - Mass Spectrometry) reference method.

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay.

Controls have to be performed :
- prior to assaying patient samples,
- at least once per day,
- after every calibration,

- and/or in accordance with laboratory and regulatory requirements.

Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PERFORMANCES

Performances were obtained on Selectra ProM, following CLSI technical recommendations, under controlled environmental conditions.

- Measuring range

a) Serum/Plasma

1.50 - 25.00 mg/dL (89 - 1487 µmol/L).

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed.

This procedure extends the measuring range up to 125.00 mg/dL (7436 µmol/L).

Do not report results outside this extended range.

b) Urine

5.0 - 250.0 mg/dL (0.30 - 14.87 mmol/L).

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 800.0 mg/dL (47.59 mmol/L).

Do not report results outside this extended range.

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

Limit of Detection (LoD) and Limit of Quantification (LoQ)

a) Serum/Plasma

LoD = 0.09 mg/dL (5 µmol/L)

LoQ = 1.00 mg/dL (59 µmol/L)

b) Urine

LoD = 0.6 mg/dL (0.04 mmol/L)

LoQ = 5.0 mg/dL (0.30 mmol/L)

- Precision

Imprecision data has been obtained on 2 Selectra ProM analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented below :

<i>a) Serum/Plasma</i>					
		Mean	Within-run	Total	
	n	mg/dL	µmol/L	CV (%)	
Level 1	80	2.41	143	0.5	2.8
Level 2	80	4.95	294	0.7	2.3
Level 3	80	6.86	408	0.7	2.2

b) Urine

		Mean	Within-run	Total	
	n	mg/dL	mmol/L	CV (%)	
Level 1	80	10.3	0.61	1.8	6.6
Level 2	80	23.9	1.42	1.1	3.8
Level 3	80	77.9	4.63	1.2	3.3

- Correlation

a) Serum/Plasma

A comparative study has been performed between URIC ACID MONO SL reagent on a Selectra ProM analyzer and a similar commercially available system on 100 human serum samples.

The sample concentrations ranged from 1.55 and 23.94 mg/dL (92 and 1424 µmol/L). The results are as follows :
Correlation coefficient : (r) = 0.999
Linear regression: y = 1.044 x - 0.04 mg/dL (2 µmol/L)

b) Urine

A comparative study has been performed between URIC ACID MONO SL reagent on a Selectra ProM analyzer and a similar commercially available system on 49 human urine samples.

The sample concentrations ranged from 5.6 and 220.2 mg/dL (0.33 and 13.10 mmol/L).

The results are as follows :
Correlation coefficient : (r) = 0.996

Linear regression: y = 1.061 x + 0.1 mg/dL (0.01 mmol/L)

- Limitations/Interferences

a) Serum/Plasma

- Studies have been performed to determine the level of interference from different compounds.

The following uric acid levels were tested: 2.52 mg/dL and 7.56 mg/dL.

No significant interference is defined by a recovery ±10% of the initial value.

Hemoglobin: No significant interference up to 50 mg/dL.

Unconjugated bilirubin: No significant interference up to 30.0 mg/dL (513 µmol/L).

Conjugated bilirubin: No significant interference up to 14.8 mg/dL (253 µmol/L).

Triglycerides: No significant interference up to 2095 mg/dL (23.7 mmol/L).

Turbidity: Interference occurs at all levels of Intralipid® Methyl-dopa: No significant interference up to 1.0 mg/dL.

Calcium dobesilate: Induces falsely low results on individuals taking calcium dobesilate.

Glucose: No significant interference up to 500 mg/dL (27.8 mmol/L).

Ascorbic acid: Significant interference on samples containing ascorbic acid.

- Do not use visibly turbid or hemolyzed samples.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.⁽⁵⁾

- Results can be falsely lowered by significant levels in the sample of NAC (N-Acetyl-Cysteine), NAPQI (metabolite of acetaminophene (paracetamol)) or metanziole.

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.⁽⁶⁻⁷⁾

b) Urine

- Studies have been performed to determine the level of interference from different compounds. The following uric acid levels were tested: 10.0 mg/dL and 75.0 mg/dL.

No significant interference is defined by a recovery ±10% of the initial value.

Hemoglobin: No significant interference up to 300 mg/dL.

Conjugated bilirubin: No significant interference up to 29.5 mg/dL (505 µmol/L).

Urea: No significant interference up to 5000 mg/dL (833 mmol/L).

Ascorbic acid: No significant interference up to 20 mg/dL.

Methyl-dopa: Induces falsely high results at therapeutic concentrations.

- Results can be falsely lowered by significant levels in the sample of NAC (N-Acetyl-Cysteine), NAPQI (metabolite of acetaminophene (paracetamol)) or metanziole.

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.⁽⁶⁻⁷⁾

- **On board stability/Calibration frequency**
On Board Stability: 28 days
Calibration frequency: 28 days
Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

These performances have been obtained using Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used.

☛The performances of applications not validated by VitalScientific are not warranted and must be defined by the user.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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