DESEMPENHO

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2. Sanhai, W.R. *et al.*, *Cardiac and Muscle Disease. Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 677 and appendix.

3. Apple, F. S. & Jaffe, A. S., *Cardiovascular Disease, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 614.

4. Boot, E., de Haan, L., *Psychopharmacology (Berl)*, (2000), 150, 347.

5. Schumann, G., *et al.*, *Clin. Chem. Lab. Med.*, (2002), 40, 635.

6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), 59, 263.

7. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).

8. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995)

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- Symbols used on our documentation are defined on ISO-15223-1 standard, except for some presented in the symbols glossary available on the VitalScientific Website. (Symbols glossary).

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
Coefficient de corrlation: (r) = 0.998

Droite de rgression : y = 1.012x + 2 U/L (

G-6-P: D-Glucose-6-Phosphate
G-6-PDH: Glucose-6-Phosphate Dehydrogenase.

The rate of increase of NADPH concentration measured at 340 nm is directly proportional to the enzymatic CK activity.

COMPOSITION

Reagent 1 : R1

Imidazole buffer, pH 6.1 (20 °C)		
D-Glucose	25	mmol/L
NADP	2.5	mmol/L
Hexokinase	≥ 6 800	U/L
Sodium azide	< 0.1	% (w/w)
Reagent 2 : R2		
Imidazole buffer, pH 8.9 (20 °C)		
Creatine phosphate	250	mmol/L
ADP	15.2	mmol/L
G-6-PDH	≥ 8 800	U/L
Sodium azide	< 0.1	% (w/w)

Also contains magnesium salts, *N*-acetyl cysteine, EDTA, diadenosine pentaphosphate and AMP for optimal performance.

MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- Normal saline solution (NaCl 9 g/L).
- Analyzers or semi-automatic analyzers.
- General Laboratory equipment (e.g. pipette).
- Do not use materials that are not required as indicated above.

PRECAUTIONS FOR USE AND WARNINGS

- The reagents R1 and R2 are classified as hazardous. (Imidazole)

Obtain Safety data sheet (SDS) before use for a proper handling.

- The reagents contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.

- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contamination.

- Do not interchange reagent vials from different kits.

STABILITY

Store at 2-8 °C and protect from light. Do not freeze. Do not use after expiration dates indicated on the vial labels.

On board stability : The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagents are ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).

- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

SAMPLES

Specimen ⁽²⁾

- Serum.
- Plasma (Lithium heparin)
- Using any other specimen type should be validated by the laboratory.

Warnings and precautions

Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

Storage and stability ⁽²⁾

- 8 hours at room temperature.
- 48 hours at 2-8°C.
- 1 month at -20°C.

REFERENCE VALUES ⁽⁹⁾

<i>Serum/plasma</i>	U/L	µkat/L		
Men	≤ 171	≤ 2.85		
Women	≤ 145	≤ 2.41		

Note : *The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

INSTALLATION AND USE

For use on **Selectra Pro analyzers :**

- Consult operator manual.

- Special Programming instructions: Programming special instructions is mandatory when some combinations of tests are performed together on the analyzer. Refer to Instructions For Use of ACID SOLUTION & SYSTEM CLEANING SOLUTION for adequate programming (See PIT-SOL).

PROCEDURE

Manual Procedure

<i>Wavelength :</i>	340 nm
<i>Optical path :</i>	1 cm
<i>Sample/ Reagent ratio :</i>	1:25
<i>Temperature:</i>	37 °C
<i>Read against distilled water.</i>	

Working reagent (4 volumes of R1 + 1 volume of R2)	1000 µL
Sample	40 µL

Mix and after 2 minutes incubation, read absorbance at 1 minute intervals during 3 minutes. Calculate the change of absorbances per minute (ΔA/min)

Automatic Procedure

These reagents may be used in several automatic analyzers. For Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

CALCULATION

Activity (U/L) = ΔA /min x 4 127

Conversion factor : U/L x 0.0167 = µkat/L

CALIBRATION

ELICAL 2 is traceable to IFCC reference method.

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay.

Controls have to be performed :

- prior to assaying patient samples,
- at least once per day,
- after every calibration,
- and/or in accordance with laboratory and regulatory requirements.

Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PERFORMANCES

Performances were obtained on Selectra ProM, following CLSI technical recommendations, under controlled environmental conditions.

- Measuring range

10 to 1714 U/L (0.17 to 28.57 µkat/L)

Samples having greater concentrations should be diluted 1:10 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 17 140 U/L (285.67 µkat/L).

Do not report results outside this extended range.

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

Limit of Detection (LoD) and Limit of Quantification (LoQ)

LoD = 1 U/L (0.02 µkat/L)

LoQ = 5 U/L (0.08 µkat/L)

- Precision

Imprecision data has been obtained on 2 Selectra Pro analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented below.

		Mean		Within-run	Total
	n	U/L	µkat/L	CV (%)	
Level 1	80	147	2.45	0.7	1.7
Level 2	80	406	6.77	1.1	2.4
Level 3	80	1154	19.23	1.1	3.9

- Correlation

A comparative study has been performed between CK NAC SL reagent on a Selectra ProM analyzer and a similar commercially available system on 100 human serum samples.

The sample concentrations ranged from 10 to 1819 U/L (0.17- 30.32 µkat/L).

The results are as follows :

Correlation coefficient : (r) = 0,998

Linear regression: y = 1.012 x + 2 U/L (0.03 µkat/L).

- Limitations/Interferences

- Hemolyzed samples should not be used since significant hemolysis may lead to falsely increased CK concentration because of adenylate kinase release.⁽¹⁾

- Studies have been performed to determine the level of interference from different compounds. The following CK levels were tested: 150 and 1200 U/L. No significant interference is defined by a recovery ±10% of the initial value.
Triglycerides: No significant interference up to 3000 mg/dL (33.9 mmol/L).
Unconjugated bilirubin: No significant interference up to 30.0 mg/dL (513 µmol/L).
Conjugated bilirubin: No significant interference up to 29.5 mg/dL (505 µmol/L).
Ascorbic acid: No significant interference up to 20.0 mg/dL.
Acetylsalicylic acid: No significant interference up to 200 mg/dL.

Acetaminophen : No significant interference up to 30 mg/dL.
- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.⁽⁶⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.^(7,8)

- **On board stability/Calibration frequency**
On Board Stability: 28 days
Calibration frequency: 28 days
Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

- **On board stability/Calibration frequency**
On Board Stability: 28 days
Calibration frequency: 28 days
Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

These performances have been obtained using Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used.
☛The performances of applications not validated by VitalScientific are not warranted and must be defined by the user.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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