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#### BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFIA/BIBLIOGRAFIA

1. Wu, A.H.B., *Clinical guide to laboratory test*, 4<sup>th</sup> Ed., (W.B. Saunders eds.), (2006), 648.
2. Sanhai W. R., *Protein Isoforms : Isoenzymes and Isoforms, Clinical Chemistry: Theory, Analysis, Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds.), (2010), 302 and appendix.
3. Panteghini, M. D., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 317.
4. Schumann, G. et al., *Clin. Chem. Lab. Med.*, (2002), **40**, 643.
5. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
7. Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AAC Press, (1997).
8. Young, D. S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AAC Press, (1995).

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- Symbols used on our documentation are defined on ISO-15223-1 standard, except for some presented in the symbols glossary available on the VitalScientific Website. (Symbols glossary).

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#### English - EN

#### INTENDED USE

LDH-L SL is an *in vitro* diagnostic reagent intended for the quantitative determination of Lactate dehydrogenase (LDH) in human serum and plasma samples on analyzers or semi-automatic analyzers.  
This *in vitro* diagnostic device is for professional use only.

#### CLINICAL SIGNIFICANCE <sup>(1-3)</sup>


Lactate dehydrogenase (LDH) can be found in nearly all cells of the body with highest activities in myocardium, liver, kidneys, skeletal muscles and erythrocytes. It is a tetramer present in the serum in the form of five main isoenzymes.  
LDH increases in case of acute myocardial infarction, hepatic disorders (viral hepatitis, cirrhosis), anemia (hemolytic, megaloblastic), some cancers and in all diseases involving tissue damage. Consequently, elevations of LDH in the serum are nonspecific.  
LDH measurement is indicated to help diagnose tissue damage in the body and sometimes to help monitor the progress of some pathologies.

#### LIMITATION OF USE

The quantitative assay of Lactate dehydrogenase (LDH) alone cannot be used to diagnose a disease or a specific pathology.  
The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

#### METHOD & PRINCIPLE <sup>(4)</sup>

IFCC method - Kinetic.

#### COMPOSITION

##### Reagent 1: R1

N-Methyl-D-Glucamine pH 9.4 (37 °C)  
Lithium L-Lactate 68 mmol/L  
Sodium azide < 0.1 % (w/w)

##### Reagent 2: R2

NAD 50 mmol/L

#### MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- Normal saline solution (NaCl 9 g/L).
- Analyzers or semi-automatic analyzers.
- General Laboratory equipment (e.g. pipette).
- Do not use materials that are not required as indicated above.

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- Sérum
- Plasma (héparine de lithium).

## PRECAUTIONS FOR USE AND WARNINGS

- Consult Safety Data Sheet (SDS) for a proper handling.

- Reagent R1 contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.

- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contamination.

- Do not interchange reagent vials from different kits.

### STABILITY

**Store at 2-8 °C and protect from light. Do not freeze.** Do not use after expiration dates indicated on the vial labels.

**On board stability :**

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

### PREPARATION

The reagent is ready to use.

### PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).

- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

### SAMPLES

#### Specimen <sup>(1,5)</sup>

- Serum

- Plasma (lithium heparin)

- Using any other specimen type should be validated by the laboratory.

**Warnings and precautions**

Samples must be free from hemolysis.<sup>(1-3)</sup>

Samples must be separated from cells and clot promptly (the presence of cells can falsely increase the result).<sup>(1-3)</sup>

Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

**Storage and stability <sup>(6)</sup>**

- 7 days at room temperature

- 4 days at 2-8°C

- 6 weeks at -20°C

In some pathologies involving LDH-4 or LDH-5 isoenzyme increase, it is preferred to store the samples at room temperature, these isoenzymes being sensitive to cold.<sup>(6)</sup>

### REFERENCE VALUES <sup>(2,3)</sup>

<i>Serum/plasma</i>	U/L	µkat/L
Adults	125 - 220	2.08 - 3.67

Higher concentrations are observed in children and new-born.

***Note :** The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

### INSTALLATION AND USE

**For use on Selectra Pro analyzers :**

- Consult operator manual.

- **Special Programming instructions: Programming special instructions is mandatory when some combinations of tests are performed together on the analyzer.** Refer to Instructions For Use of ACID SOLUTION & SYSTEM CLEANING SOLUTION for adequate programming (See PIT-SOL).

### PROCEDURE

<b>Manual Procedure</b>		
<i>Wavelength<span> </span>:</i>	340 nm	
<i>Optical path<span> </span>:</i>	1 cm	
<i>Sample/ Reagent ratio<span> </span>:</i>	1:36	
<i>Temperature<span> </span>:</i>	37 °C	
<i>Read against reagent blank.</i>		
	BLANK	SAMPLE
<b>Reagent R1</b>	1000 µL	1000 µL
<b>Distilled water</b>	35 µL	-
<b>Sample</b>	-	35 µL

*Mix, wait 3 minute and add:*

<b>Reagent R2</b>	250 µL	250 µL
-------------------	--------	--------

*Mix and after a 100 second incubation, read absorbance at 1 minute intervals during 3 minutes. Calculate the change of absorbance per minute (ΔA/min).*

#### Automatic Procedure

These reagents may be used on several automatic analyzers. For Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

### CALCULATION

Activity (U/L) = ΔA/min x 5828

**Conversion factor :** U/L x 0.0167 = µkat/L

### CALIBRATION

ELICAL 2 is traceable to IFCC reference method.

**Calibration frequency :** The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

### QUALITY CONTROL

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay.

Controls have to be performed :

- prior to assaying patient samples,

- at least once per day,

- after every calibration,

- and/or in accordance with laboratory and regulatory requirements.

Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

### WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

### PERFORMANCES

Performances were obtained on Selectra ProM, following CLSI technical recommendations, under controlled environmental conditions.

### - Measuring range

50 - 800 U/L (0.83 - 13.33 µkat/L)

Samples having greater concentrations should be diluted 1:10 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 8000 U/L (133.33 µkat/L).

Do not report results outside this extended range.

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

#### Limit of Detection (LoD) and Limit of Quantification (LoQ)

LoD = 4 U/L (0.07 µkat/L)

LoQ = 10 U/L (0.17 µkat/L)

#### - Precision

Imprecision data has been obtained on 2 Selectra ProM analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented below.

		Mean	Within-run	Total	
	n	U/L	µkat/L	CV (%)	
Level 1	80	168	2.80	0.7	4.1
Level 2	80	309	5.15	0.7	2.8
Level 3	80	712	11.87	0.5	3.0

#### - Correlation

A comparative study has been performed between LDH-L SL reagent on a Selectra ProM analyzer and a similar commercially available system on 99 human serum samples.

The sample concentrations ranged from 45 to 780 U/L (0.75 - 13.00 µkat/L).

The results are as follows :

Correlation coefficient : (r) = 0.997

Linear regression: y = 1.010 x + 3 U/L (0.05 µkat/L)

#### - Limitations/Interferences

- Hemolyzed samples should not be used since significant hemolysis may lead to falsely increased results due to high levels of LDH originating from erythrocytes.<sup>(1-3)</sup>

- Studies have been performed to determine the level of interference from different compounds.

The following Lactate dehydrogenase (LDH) levels were tested: 200 U/L and 700 U/L

No significant interference is defined by a recovery ±10% of the initial value.

**Unconjugated bilirubin:** No significant interference up to 30.0 mg/dL (513 µmol/L).

**Conjugated bilirubin:** No significant interference up to 29.5 mg/dL (505 µmol/L).

**Triglycerides:** No significant interference up to 3146 mg/dL (35.6 mmol/L).

**Ascorbic acid:** No significant interference up to 20.0 mg/dL.

**Acetylsalicylic Acid:** No significant interference up to 200 mg/dL.

**Acetaminophen:** No significant interference up to 30 mg/dL.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström macroglobulinemia) can cause unreliable results.<sup>(6)</sup>

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.<sup>(7-8)</sup>

#### - On board stability/Calibration frequency

**On Board Stability:** 28 days

**Calibration frequency:** 6 days

Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

*These performances have been obtained using Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used.*

*☛The performances of applications not validated by VitalScientific are not warranted and must be defined by the user.*

### DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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