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4. Daly J.A. & Erlingshausen G., *Clin. Chem.*, (1972), 18, 263.

5. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature. Acta Clin Belg.*, (2004), 59, 263.

6. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACCPress, (1997).

7. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACCPress, (1995).

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- Symbols used on our documentation are defined on ISO-15223-1 standard, except for some presented in the symbols glossary available on the VitalScientific Website. (Symbols glossary).

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- **Special Programming instructions** : see § INSTALLATION AND USE

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CLÍNICAL SIGNIFICANCE ⁽¹⁻³⁾

Phosphorus is present throughout the human body as organic or inorganic phosphate: 80-85% in bones (hydroxyapatite), 15% in soft tissues (mainly in cells, integrated in nucleic acids, phospholipids, phospho-proteins, ATP, phosphorylated carbohydrates), less than 0.1% in extracellular fluids. In blood, phosphorus is mainly present as inorganic phosphate (Pi).

Levels of Pi are closely linked to that of Calcium. An elevation of serum Pi can occur in renal failure, hyperparathyroidism or pseudo-hyperparathyroidism, excessive phosphate intake, leakage of intra-cellular phosphorus (e.g. untreated diabetic ketoacidosis) and vitamin D intoxication. A decrease of serum Pi can occur in hyperparathyroidism, renal tubule defects (Fanconi's syndrome), vitamin D deficiency (rickets, osteomalacia) and diabetic ketoacidosis induced by insulin.

Measuring Pi in urine can be recommended as an additional test to serum Pi.

In practice, measurement of Pi in blood or urine is indicated to help diagnose pathologies involving abnormal serum calcium, for the monitoring of kidney pathologies, for the follow-up of diabetic patients or to assess gastrointestinal disorders.

LIMITATION OF USE

The quantitative assay of inorganic phosphorus alone cannot be used to diagnose a disease or a specific pathology.

The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

METHOD & PRINCIPLE ⁽⁴⁾

Phosphomolybdate - End point

Inorganic Phosphate + H₂SO₄ + Ammonium Molybdate → Phosphomolybdate Complex

COMPOSITION

Reagent: R

Sulfuric acid 210 mmol/L

Ammonium molybdate 650 µmol/L

Surfactant

Standard: Std

Phosphorus 4.7 mg/dL

1.52 mmol/L

MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2

- CONT-0060 ELITROL I

- CONT-0160 ELITROL II

- Normal saline solution (NaCl 9 g/L).

- Analyzers or semi-automatic analyzers.

- General Laboratory equipment (e.g. pipette).

- Do not use materials that are not required as indicated above.

PRECAUTIONS FOR USE AND WARNINGS

- The reagent R is classified as hazardous.

- Consult Safety Data Sheet (SDS) for a proper handling.

- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contamination.

STABILITY

Store at 2-25 °C and protect from light. Do not freeze.

Do not use after expiration dates indicated on the vial labels.

The standard should be immediately and tightly capped to prevent contamination and evaporation.

On board stability :

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagent and standard are ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).

- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

SAMPLES

Specimen ^(1,2)

Serum

- Plasma (lithium heparin)

- Urine

- Using any other specimen type should be validated by the laboratory.

Warnings and precautions

- Collection of blood samples must be performed on fasting patients. ⁽³⁾

- Separate blood samples from erythrocytes within 1 hour after collection. ⁽³⁾

- Blood samples must be free from hemolysis. ⁽¹⁾

- Urines must be acidified with HCl after collection (pH<3). ⁽³⁾

- Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.

Storage and stability ⁽²⁾

Serum/ plasma

- 4 days at room temperature

- 7 days at 2-8°C

- Several months at -20°C

Urine

- 6 months

REFERENCE VALUES ⁽¹⁾

<i>Serum/plasma</i>	mg/dL	mmol/L	
Adults	2.5 - 4.5	0.81 - 1.45	
Children	4.0 - 7.0	1.29 - 2.26	

<i>Urine (24h collection)</i>	g/24h	mmol/24h	
	0.4 - 1.3	12.9 - 42.0	
	For a urinary volume of 1.5 L/24h		
	mg/dL	mmol/L	
	26.7 - 86.7	8.6 - 28.0	

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

INSTALLATION AND USE

For use on *Selectra Pro analyzers :*

- Consult operator manual.

- **Special Programming instructions: Programming special instructions is mandatory when some combinations of tests are performed together on the analyzer.** Refer to Instructions For Use of ACID SOLUTION & SYSTEM CLEANING SOLUTION for adequate programming (See PIT-SOL).

PROCEDURE

Manual Procedure

Wavelength : 340 nm

Optical path : 1 cm

Sample/ Reagent ratio : 1:100

Temperature : 37 °C

Urine samples must be diluted 1:10 with NaCl 9 g/L solution before measurement

Read against reagent blank.

	CALIBRATION	TEST
Reagent R	1 000 µL	1 000 µL
Standard/ Calibrator	10 µL	-
Sample	-	10 µL

Mix and read the absorbances (A) after an incubation of 5 minutes.

Automatic Procedure

These reagents may be used on several automatic analyzers. For Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

Urine samples must be diluted 1:10 with NaCl 9 g/L solution before measurement. For users of Selectra TouchPro software, urine dilution is performed automatically.

CALCULATION

^A Sample x n n = Standard/Calibrator concentration

^A Standard/ Calibrator

For the calculation of phosphorus concentration in urine, multiply the result by the dilution factor (10). For users of Selectra TouchPro software, the results take the dilution factor into account.

Conversion factor : mg/dL x 0.3229 = mmol/L

CALIBRATION

ELICAL 2 and Phosphorus standard are traceable to SRM 956d reference material.

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay.

Controls have to be performed :

- prior to assaying patient samples,

- at least once per day,

- after every calibration,

- and/or in accordance with laboratory and regulatory requirements.

Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PERFORMANCES

Performances were obtained on Selectra ProM, following CLSI technical recommendations, under controlled environmental conditions.

- Measuring range

a) Serum/Plasma

2.00 - 20.00 mg/dL (0.65 - 6.46 mmol/L).

Do not report results outside the measuring range.

b) Urine

10.0 - 200.0 mg/dL (3.2 - 64.6 mmol/L).

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 1 000.0 mg/dL (322.9 mmol/L).

Do not report results outside this extended range.

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

Limit of Detection (LOD) and Limit of Quantification (LoQ)

a) Serum/Plasma

LoD = 0.12 mg/dL (0.04 mmol/L)

LoQ = 1.00 mg/dL (0.32 mmol/L)

b) Urine

LoD = 0.2 mg/dL (0.1 mmol/L)

LoQ = 10.0 mg/dL (3.2 mmol/L)

- Precision

Imprecision data has been obtained on 2 Selectra ProM analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented below :

<i>a) Serum/Plasma</i>					
		Mean	Within-run	Total	
	n	mg/dL	mmol/L	CV (%)	
Level 1	80	2.39	0.77	0.6	2.1
Level 2	80	4.77	1.54	0.4	2.3
Level 3	80	9.02	2.91	0.7	2.3

b) Urine

		Mean	Within-run	Total	
	n	mg/dL	mmol/L	CV (%)	
Level 1	80	25.2	8.1	1.0	2.0
Level 2	80	72.5	23.4	0.6	1.0
Level 3	80	144.2	46.6	0.5	1.2

- Correlation

a) Serum/Plasma

A comparative study has been performed between PHOSPHORUS reagent on a Selectra ProM analyzer and a similar commercially available system on 100 human serum samples. The sample concentrations ranged from 2.01 to 18.82 mg/dL (0.65 - 6.08 mmol/L).

The results are as follows :

Correlation coefficient : (r) = 0.999

Linear regression: y = 1.004 x + 0.05 mg/dL (0.02 mmol/L)

b) Urine

A comparative study has been performed between PHOSPHORUS reagent on a Selectra ProM analyzer and a similar commercially available system on 52 human urine samples.

The sample concentrations ranged from 10.2 to 182.3 mg/dL (3.3 - 58.9 mmol/L).

The results are as follows :

Correlation coefficient : (r) = 0.999

Linear regression: y = 0.883 x - 0.3 mg/dL (0.1 mmol/L)

- Limitations/Interferences

a) Serum/Plasma

- Hemolyzed samples should not be used since significant hemolysis may falsely increase the result due to the high level of phosphates in erythrocytes. ⁽¹⁾

- Studies have been performed to determine the level of interference from different compounds. The following inorganic phosphorus levels were tested: 2.17 mg/dL and 4.65 mg/dL.

No significant interference is defined by a recovery ≤±10% of the initial value.

Unconjugated bilirubin : No significant interference up to 15.0 mg/dL (257 µmol/L).

Conjugated bilirubin : No significant interference up to 3.0 mg/dL (51 µmol/L).

Glucose : No significant interference up to 528 mg/dL (29.3 mmol/L).

Triglycerides : No significant interference up to 729 mg/dL (8.2 mmol/L).

Hemoglobin : No significant interference up to 50 mg/dL.

Ascorbic acid : No significant interference up to 20.0 mg/dL.

Methyl-Dopa : No significant interference up to 1.0 mg/dL.

Do not use visibly turbid, icteric or hemolyzed samples.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results. ⁽²⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young. ⁽⁶⁻⁷⁾

b) Urine

- Studies have been performed to determine the level of interference from different compounds.

The following inorganic phosphorus levels were tested: 25.0 mg/dL and 150.0 mg/dL.

No significant interference is defined by a recovery ≤±10% of the initial value.

Conjugated bilirubin : No significant interference up to 22.1 mg/dL (378 µmol/L).

Hemoglobin : No significant interference up to 500 mg/dL.

Uric acid : No significant interference up to 120 mg/dL (7.1 mmol/L).

Urea : No significant interference up to 6 000 mg/dL (999 mmol/L).

Ascorbic acid : No significant interference up to 20.0 mg/dL.

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young. ⁽⁶⁻⁷⁾

- **On board stability/Calibration frequency**

On Board Stability: 28 days

Calibration frequency: 28 days

Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

These performances have been obtained using Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used.

☛The performances of applications not validated by VitalScientific are not warranted and must be defined by the user.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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