

▶ Scalable Solutions

VitalScientific offers a comprehensive range of bench-top chemistry systems for the physician office laboratory and small hospitals. Our bench-top chemistry systems are designed to meet the needs of laboratories running from 15,000 to 200,000 tests per year. All systems are supported by technical service and remote diagnostics.

Leading Family of Bench-top Chemistry Systems



INSTRUMENT	PATIENTS PER DAY	TESTS PER HOUR	ANNUAL TEST VOLUME
Selectra ProS	5 - 15	140	15,000 - 60,000
Selectra ProM	10 - 40	180	40,000 - 100,000
ENVOY 500+	20 - 80	490	50,000 - 200,000



Selectra ProS™

Fully-featured, compact system

- 5-15 patients per day
- 30 refrigerated reagent positions
- 25 position sample rotor



Selectra ProM™

Fully-featured, mid-volume system

- 10-40 patients per day
- 32 refrigerated reagent positions
- 50 position sample rotor



ENVOY® 500+

Nothing faster in bench-top chemistry

- 20-80 patients per day
- 40 refrigerated reagent positions
- 52 position sample wheel
- Clot detection sensor
- Reusable glass cuvettes

Clinical Chemistry

Standard features for all analyzers

- Touch screen user interface
- Ready-to-use, liquid-stable reagents
- Positive sample identification, host query¹
- Barcoded reagents
- Liquid level sensing
- Programmable result checks, auto-dilutions
- Remote diagnostics

► Comprehensive Chemistry Systems Service and Support

- Technical Support, Monday – Friday, 8AM – 7PM EST.
- Remote Diagnostics and on-site service
- Full one-year warranty; post-warranty options

► VitalScientific

VitalScientific is a global leader in benchtop Clinical Chemistry, with over 50,000 analytical systems delivered worldwide. Our Scalable Systems prioritize reliability and analytical quality, empowering better medical decisions closer to the patient. With 20+ years servicing Physician Office Laboratories (POL), clinics, and small reference labs in the USA, we're renowned for excellent product quality and user-friendliness. ENVOY® 500/500+ and Selectra Pro Series offer high-quality solutions, including dedicated barcoded reagents, calibrators, controls, and consumables, supporting laboratories in delivering top-tier healthcare. We facilitate the laboratories in their key role of providing first-class healthcare.

VitalScientific
info@vitalscientific.com

©2024, VitalScientific, The Netherlands



VitalScientific
Diagnostics for Life

www.vitalscientific.com